SINGLE-USE BIPOLAR FORCEPS

- DE** GEBRAUCHSANWEISUNG
- EN** INSTRUCTIONS FOR USE
- US** INSTRUCTIONS FOR USE
- PT** INSTRUÇÕES DE USO
- ES** INSTRUCCIONES DE USO
- IT** ISTRUZIONI PER L'USO

DEUTSCH

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
ENGLISH

INSTRUCTIONS FOR USE SINGLE-USE BIPOLAR FORCEPS

REF Product group: Single-Use HF Forceps

adeor medical AG
Martinshof 5
83626 Valley
Germany

Attention: The single-use instruments are supplied sterile.

Rx Only Warning: Legal regulations stipulate that the single-use instruments may only be acquired by or on behalf of a physician.

Warning: Please read the information in this leaflet carefully. Improper use and care or use for unintended purposes can lead to premature wear or risks for patients and users.

Intended use: Single-use bipolar forceps made by adeor medical AG are intended for use in all areas of open surgery. They are designed to grasp, manipulate and coagulate tissue. They are connected through a suitable bipolar cable with the bipolar output of a HF generator and may be used only with bipolar coagulation current. adeor bipolar forceps must be operated with the following parameters: Frequency range between 300 kHz and 1,000 kHz; maximum generator operating voltage 600 Vp. The surgeon must determine the appropriate type of treatment for each individual case in cooperation with the internist and the anaesthetist. Typical HF-surgical interventions include

- Dermatology: hair and wart removal
- ENT: polyp removal, bipolar hemostasis
- Gynecology: electro-conization of the portio, sterilization by bipolar coagulation of the fallopian tube
- Urology: transurethral resection of the prostate, bladder tumor or bladder wall papillomas; cutting and coagulation in transabdominal surgery
- Laparoscopy and cystoscopy: Removal of polyps
- General surgery: primary sections, separating tissue, plastic surgery, hemostasis, coagulation
- Neurosurgery: Bipolar coagulations

Contraindications: The single-use bipolar forceps should not be used if the physician determines that the risk to the patient outweighs the benefits. Single-Use bipolar forceps made by adeor medical AG should not be used for tubal sterilization or tube coagulation for sterilization. Do not use on patients with cardiac pacemakers or other active implants.

Attention: Single-use instruments for electrosurgery may only be used by persons who are specially trained or instructed.

Safety instructions:

- Checking, using and handling electrosurgical instruments is the responsibility of the user.
- Before each use, the single-use instrument must be inspected for wear and visible damage, such as cracks, fractures or damage to the insulation, in particular areas such as the blade, tip, groove, connector, lock and moving parts. Check the insulation thoroughly.
- Do not use damaged single-use instruments.
- Do not use the single-use instruments in the presence of flammable or explosive substances.

- Do not place the single-use instrument on top of the patient.
- Perform coagulation only if the contact surfaces are in the visible area. Do not touch any other metal instruments during the procedure.

Warning: The parts mentioned in this instruction for use are **single-use only**. They cannot be resterilized. Resterilization causes malfunction to parts of the product which infects the function. If the single-use forceps are damaged, especially at the insulation or the connector, dispose the product.

Single-Use (Non-stick) Bipolar forceps: The polished tips are made of silver (with gold plating) and may tarnish which **does not** affect the function of the single-use instruments.

Functional testing: Single-use instruments must be checked for functionality before each use. Visually inspect the instruments for cleanliness and carry out a function test. If necessary, repeat the reprocessing process until the instrument looks clean. Surface damage, such as scratches, cracks, nicks, notches etc., or bent parts mean that the instruments may not be used.

Warning: Malfunction, instrument failure, damaged insulation

Action required by user: Perform a visual inspection/functional testing. It is the responsibility of the user to ensure that the instruments are used in line with the intended purpose.

Warning: Electrical loading capacity of the electrodes exceeded

Action required by user: Follow manufacturer's instructions for electrical loading; see "Intended use".

Handling: Surgical instruments should always be treated with the utmost care. This is particular pointed on tips and other sensitive areas.

Generator connections: Single-use forceps made by adeor medical AG are approved for use with the following generators:

ERBE Elektromedizin GmbH		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
VIO 300 D	10140-100	≤ 200 Watt
VIO 200 D	10140-200	≤ 120 Watt
VIO 300 S	10140-300	≤ 200 Watt
VIO 200 S	10140-400	≤ 120 Watt
VIO 100 C	10140-500	≤ 80 Watt
VIO 50 C	10140-550	≤ 50 Watt
ICC 350	ICC 350	≤ 150 Watt
ICC 300	ICC 300	≤ 150 Watt
ICC 200	ICC 200	≤ 150 Watt
ICC 80	ICC 80	≤ 50 Watt
ICC 50	ICC 50	≤ 50 Watt
ACC 451	ICC 451	≤ 120 Watt
ACC 450	ICC 450	≤ 120 Watt
ACC 430	ICC 43	≤ 120 Watt

SUTTER Medizintechnik GmbH		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
BM-780 II	360080-01	≤ 70 Watt

KARL STORZ GmbH & Co. KG		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
AUTOCON® II 200	205322 20	≤ 120 Watt
AUTOCON® II 400	205352 20	≤ 120 Watt

COVIDIEN		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
Force FX™	Force FX™	≤ 70 Watt
ForceTriad™	ForceTriad™	≤ 95 Watt
Force EZ™	Force EZ™	≤ 70 Watt
SurgiStat™	SurgiStat™	≤ 30 Watt

BOWA-electronic GmbH & Co. KG		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
ARC 400	900-400	≤ 350 Watt
ARC 350	900-351	≤ 350 Watt
ARC 303	900-303	≤ 120 Watt
ARC 250	900-250	≤ 120 Watt
ARC 100	900-100	≤ 100 Watt

Olympus Surgical Technologies Europe		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
ESG-100, 220-240 V	WB991036	≤ 120 Watt
ESG-100, 100-120 V	WB991046	≤ 120 Watt
ESG-400	WB91051W	≤ 320 Watt

Söring GmbH		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
MBC 600	MBC 600	≤ 100 Watt
MBC 601	MBC 601	≤ 100 Watt
UAM	UAM	
MBC 200	MBC 200	≤ 140 Watt
BCC 140	BCC 140	≤ 140 Watt

LAMIDEY NOURY MEDICAL		
Typ	ITEM NO.	MAX. BIPOLAR COAG. POWER
Optima 4	Optima 4	≤ 370 Watt
Optima 3	Optima 3	≤ 370 Watt
Optima 2	Optima 2	≤ 200 Watt
MC 2	MC 2	≤ 70 Watt
MC 3	MC 3	≤ 100 Watt
MC 4	MC 4	≤ 100 Watt

Integra LifeSciences Corporation		
Typ	ITEM NO.	MAX.

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Gebrüder Martin GmbH & Co. KG
Aesculap AG / B. Braun Melsungen AG