SINGLE-USE BIPOLAR FORCEPS

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EN INSTRUCTIONS FOR USE

US INSTRUCTIONS FOR USE

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DEUTSCH

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
ENGLISH

INSTRUCTIONS FOR USE SINGLE-USE BIPOLAR FORCEPS

REF	Product group: Single-Use HF Forceps
------------	---

Reference: BIPAAABBSU
AAA will be replaced by a 3-digit number.
BB will be replaced by EU or US.
EU represents the Euro-pean and US the American connector type.
SU means Single-Use.

adeor medical AG Martinshof 5 83626 Valley Germany

Attention: The single-use instruments are supplied sterile.

Rx Only Warning: Legal regulations stipulate that the single-use instruments may only be acquired by or on behalf of a physician.

Warning:

Please read the information in this leaflet carefully. Improper use and care or use for unintended purposes can lead to premature wear or risks for patients and users.

Intended use:

Single-use bipolar forceps made by adeor medical AG are intended for use in all areas of open surgery. They are de-signed to grasp, manipulate and coagulate tissue. They are connected through a suitable bipolar cable with the bipolar output of a HF generator and may be used only with bipolar coagulation current. adeor bipolar forceps must be opera-ted with the following parameters:
Frequency range between 300 kHz and 1,000 kHz; maxi-mum generator operating voltage 600 Vp.
The surgeon must determine the appropriate type of treat-ment for each individual case in cooperation with the internist and the anaesthetist. Typical HF surgical interventions include

- Dermatology: hair and wart removal
- ENT: polyp removal, bipolar hemostasis
- Gynecology: electro-conization of the portio, sterilization by bipolar coagulation of the fallopian tube
- Urology: transurethral resection of the prostate, bladder tumour or bladder wall papillomas; cutting and coagulation in transabdominal surgery
- Laparoscopy and cystoscopy: Removal of polyps
- General surgery: primary sections, separating tissue, plastic surgery, hemostasis, coagulation
- Neurosurgery: Bipolar coagulations

Contraindications:

The single-use bipolar forceps should not be used if the physician determines that the risk to the patient outweighs the benefits.

Single-Use bipolar forceps made by adeor medical AG should not be used for tubal sterilization or tube coagulation for sterilization.

Do not use on patients with cardiac pacemakers or other active implants.

Attention:

Single-use instruments for electrosurgery may only be used by persons who are specially trained or instructed.

Safety instructions:

- Checking, using and handling electrosurgical instruments is the responsibility of the user.
- Before each use, the single-use instrument must be inspected for wear and visible damage, such as cracks, fractures or damage to the insulation, in particular areas such as the blade, tip, groove, connector, lock and moving parts. Check the insulation thoroughly.
- Do not use damaged single-use instruments.
- Do not use the single-use instruments in the presence of flammable or explosive substances.

- Do not place the single-use instrument on top of the patient.
- Perform coagulation only if the contact surfaces are in the visible area. Do not touch any other metal instru-ments during the procedure.

Warning:

The parts mentioned in this instruction for use are **single-use only**. They cannot be resterilized.

Resterilization causes malfunction to parts of the product which infects the function.

If the single-use forceps are damaged, especially at the insulation or the connector, dispose the product.

Single-Use (Non-stick) Bipolar forceps: The polished tips are made of stainless steel (silver) and may tarnish which **does not** affect the function of the single-use instruments.

Functional testing:

Single-use instruments must be checked for functionality before each use. Visually inspect the instruments for clean-liness and carry out a function test. If necessary, repeat the reprocessing process until the instrument looks clean. Surface damage, such as scratches, cracks, nicks, notches etc., or bent parts mean that the instruments may not be used.

Warning:

Malfunction, instrument failure, damaged insulation

Action required by user: Perform a visual inspection/functional testing. It is the re-sponsibility of the user to ensure that the instruments are used in line with the intended purpose.

Warning:

Electrical loading capacity of the electrodes exceeded

Action required by user: Follow manufacturer's instructions for electrical loading; see "Intended use".

Handling:

Surgical instruments should always be treated with the utmost care. This is particular pointed on tips and other sensitive areas.

Generator connections:

Single-use forceps made by adeor medical AG are appro-ved for use with the following generators:

ERBE Elektromedizin GmbH		
TYPE	ITEM NO.	MAX. BIPOLAR COAG. POWER
VIO 300 D	10140-100	≤ 200 Watt
VIO 200 D	10140-200	≤ 120 Watt
VIO 300 S	10140-300	≤ 200 Watt
VIO 200 S	10140-400	≤ 120 Watt
VIO 100 C	10140-500	≤ 80 Watt
VIO 50 C	10140-550	≤ 50 Watt
ICC 350	ICC 350	≤ 150 Watt
ICC 300	ICC 300	≤ 150 Watt
ICC 200	ICC 200	≤ 150 Watt
ICC 80	ICC 80	≤ 50 Watt
ICC 50	ICC 50	≤ 50 Watt
ACC 451	ICC 451	≤ 120 Watt
ACC 450	ICC 450	≤ 120 Watt
ACC 430	ICC 43	≤ 120 Watt

SUTTER Medizintechnik GmbH		
TYPE	ITEM NO.	MAX. BIPOLAR COAG. POWER
BM-780 II	360080-01	≤ 70 Watt

KARL STORZ GmbH & Co. KG		
TYPE	ITEM NO.	MAX. BIPOLAR COAG. POWER
AUTOCON® II 200	205322 20	≤ 120 Watt
AUTOCON® II 400	205352 20	≤ 120 Watt

COVIDIEN		
TYPE	ITEM NO.	MAX. BIPOLAR COAG. POWER
Force FX™	Force FX™	≤ 70 Watt
ForceTriad™	ForceTriad™	≤ 95 Watt
Force EZ™	Force EZ™	≤ 70 Watt
SurgiStat™	SurgiStat™	≤ 30 Watt

BOWA-electronic GmbH & Co. KG		
TYPE	ITEM NO.	MAX. BIPOLAR COAG. POWER
ARC 400	900-400	≤ 350 Watt
ARC 350	900-351	≤ 350 Watt
ARC 303	900-303	≤ 120 Watt
ARC 250	900-250	≤ 120 Watt
ARC 100	900-100	≤ 100 Watt

US

INSTRUCTIONS FOR USE
SINGLE-USE BIPOLAR FORCEPS

REF	Product group: Single-Use HF Forceps
------------	---

Reference: BIPAAABBSU
AAA will be replaced by a 3-digit number.
BB will be replaced by EU or US.
EU represents the European and US the American connector type.
SU means Single-Use.

 adeor medical AG
Martinhof 5
83626 Valley
Germany

Attention:

The single-use instruments are supplied sterile.

Warning:

Legal regulations stipulate that the single-use instruments may only be acquired by or on behalf of a physician.

Warning:

Please read the information in this leaflet carefully. Improper use and care or use for unintended purposes can lead to premature wear or risks for patients and users.

Intended use:

The adeor medical Single-Use Non-stick Bipolar Forceps are intended for use by a physician familiar with electro-surgery in bipolar coagulation for general open surgery where coagulation of soft tissue is needed.
adeor bipolar forceps must be operated with the following parameters:
Frequency range between 300 kHz and 1,000 kHz;
maximum generator operating voltage 600 Vp.
The adeor medical Single-Use Non-stick Bipolar Forceps have not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures and should not be used for these procedures.

The types of surgery intended include:

- ENT
- Gynecology (except for use in female sterilization)
- Urology
- General surgery
- Neurosurgery
- Laryngeal Surgery
- Orthopedic Surgery
- Thoracic Surgery

Contraindications:

The single-use bipolar forceps should not be used if the physician determines that the risk to the patient outweighs the benefits.
Single-Use bipolar forceps made by adeor medical AG should not be used for tubal sterilization or tube coagulation for sterilization.

Do not use on patients with cardiac pacemakers or other active implants.

Attention:

Single-use instruments for electrosurgery may only be used by persons who are specially trained or instructed.

Safety instructions:

- Checking, using and handling electrosurgical instruments is the responsibility of the user.
- Before each use, the single-use instrument must be inspected for wear and visible damage, such as cracks, fractures or damage to the insulation, in particular areas such as the blade, tip, groove, connector, lock and moving parts. Check the insulation thoroughly.
- Do not use damaged single-use instruments.
- Do not use the single-use instruments in the presence of flammable or explosive substances.
- Do not place the single-use instrument on top of the patient.
- Perform coagulation only if the contact surfaces are in the visible area. Do not touch any other metal instruments during the procedure.

Warning:

The parts mentioned in this instruction for use are single-use only. They cannot be resterilized.

Resterilization causes malfunction to parts of the product which infects the function.

If the single-use forceps are damaged, especially at the insulation or the connector, dispose the product.

Single-Use (Non-stick) Bipolar forceps:

The polished tips are made of stainless steel (silver) and may tarnish which does not affect the function of the single-use instruments.

Functional testing:

Single-use instruments must be checked for functionality before each use. Visually inspect the instruments for cleanliness and carry out a function test. If necessary, repeat the reprocessing process until the instrument looks clean.

Surface damage, such as scratches, cracks, nicks, notches etc., or bent parts mean that the instruments may not be used.

Warning:

Malfunction, instrument failure, damaged insulation

Action required by user:

Perform a visual inspection/functional testing. It is the responsibility of the user to ensure that the instruments are used in line with the intended purpose.

Warning:

Electrical loading capacity of the electrodes exceeded

Action required by user:

Follow manufacturer's instructions for electrical loading; see "Intended use".

Handling:

Surgical instruments should always be treated with the utmost care. This is particular pointed on tips and other sensitive areas.

Suitable connection cables:

adeor single-use bipolar forceps can be used with bipolar connection cables that are suitable for use with generators that meet the required operating parameters.

Cable connection, forceps

European flat plug

2-pin plug (US male connector)

These cables may be connected only to the bipolar output of electrosurgical devices. Please follow the generator instructions on connecting these cables before using the single-use bipolar forceps.

HF generators pose a number of risks, such as incorrect operation, unintended high-frequency burns, ignition of flammable liquids and gases (risk of explosion).

Packing/Storage:

Some parts of these single-use instruments can be easily damaged. The sterilized single-use instruments should be stored in a dry, clean and dust-free environment at moderate temperatures of 5°C to 40°C.

Warranty:

adeor medical AG provides customers with tested, fault-free products. All our products are designed and manufactured to meet highest quality standards. However, should problems occur, please contact our customer service team. Liability is excluded for products that were modified from the original product, used for other than the intended purposes, handled or used improperly.

adeor medical AG does not accept liability for incidental or consequential damages. This warranty is voided if repairs are carried out by companies that are not authorized for repair by adeor medical AG.

Proper disposal of bipolar forceps:

If the single-use bipolar forceps can no longer be used, e.g., because of damage, they must be disposed properly.

(This section is a translation of the German text. It may not be accurate.)

Symbols and their meaning:

	Follow operating instructions
	Protect the product against moisture and condensation
	Protect the product against direct sunlight


	Sterilized by radiation
---	-------------------------

	Caution: Please follow instructions
---	-------------------------------------

	Manufacturer / Date of manufacture
---	------------------------------------


	Order No.
--	-----------

	Do not reuse/Single-use only
---	------------------------------

	Lot number
---	------------

	Use by – symbol is accompanied by a date
---	--

	Prescription device only (USA)
---	--------------------------------

	CE including Notified Body code
---	---------------------------------

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Söring GmbH
Integra LifeSciences Corporation
Gebrüder Martin GmbH & Co. KG
Aesculap AG / B. Braun Melsungen AG