CRANIOBRIDGE

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INSTRUCTIONS FOR USE CRANIOBRIDGE™ Cranial Fixation

General

Please read the following information carefully. Improper handling, maintenance, and inappropriate use can cause premature wear or risks for patients and user. adeor medical AG products may only be used with original adeor medical AG equipment.

adeor medical AG implants may be used only by qualified and certified medical staff. In case of technical issue, it is to ensure that an alternative system is available.

The products are medical implantable devices and are designated to be used on human beings. The products are not designated to be used on the central nervous system.

To guarantee the safe application of the products only adeor medical AG original equipment is allowed to be used with the product.

Intended Use

The implants consist of firm material that is attached with screws to fractured craniofacial bone (including the maxilla and/or mandible), to bridge and stabilize the fracture gap, and shield the site from stress as the bone heals; the material cannot be chemically degraded or absorbed by natural body processes. The device may also be used for craniofacial reconstructive procedures (e.g., to treat congenital malformations), endobrow fixation, or craniotomy flap fixation.

The implants, as well as the accessories, must never be used outside of this intended purpose.

Indications

- Trauma repair and reconstruction of the craniomaxillofacial skeleton
- Craniofacial reconstructive procedures
- Endobrow fixation
- Craniotomy flap fixation
- Orbital floor fractures
- Medial orbital wall fractures

Contraindications

The contraindications for the use of the implants are depending on the single disease pattern and the intended surgical procedure.

General Contraindications:

- Health condition with inhibiting healing process or insufficient implant support (e.g. interference of the blood supply; Insufficient bone quality or quantity)
- Local bone tumors
- Fever or leucocytosis
- The mental conditions which make a participation in the rehabilitation program impossible (Parkinson disease, alcoholism, drug-usage, etc.)
- Systemic diseases and metabolic disorders; infections and falls; obesity
- Pregnancy
- Severe anatomical deformities (e.g., due to birth defects)
- Need to combine metals, components or different systems
- Insufficient bone substance(s) or tissue inclusion
- Disruption of anatomical structures or limitation of physiological performance due to application
- Absence of patient compliance with postoperative instructions

Relative Contraindications:

- Rapidly progressive joint disease, bone resorption, osteopenia, osteomalacia, and/or osteoporosis

Specific Contraindications:

- Allergy to a material component
- Infections
- Acute and chronic infections; muscle, nerve or vascular diseases

Material

The implants are made of titanium according to

- ISO 5832-2 Pure titanium (unalloyed),
 - ISO 5832-3 Ti6Al4V alloy,
 - ISO 5852-11 Ti6Al7Ni alloy,
- are biocompatible and prevent so-called chromium-nickel allergies in their properties.

The surface of these implants is chemically passive and non-magnetic.

These implants can be combined with standardized materials, the composition of which is within the standardized directional analysis according to ISO 5832-2 / ISO 5832-3 / ISO 5832-11 and the required specifications.

MRI Note

adeor medical AG implants are made of titanium and titanium alloys. This type of non-ferromagnetic material poses no risk to the patient during magnetic resonance imaging (MRI) examinations.

Various studies and test reports have shown that these materials can be safely used in an MRI environment according to ISO 5832-2 and ISO 5832-3 and can be considered „MRI Conditional“ according to ASTM F2503. MRI Conditional is defined as an item that has been shown to pose no known

hazards in a specific MRI environment with specified conditions of use (ASTM F2503, 3.1.9).

Precautions:

- The above statement is based on non-clinical testing. The actual patient temperature rise depends on a variety of factors beyond SAR and timing of radiofrequency application. Therefore, it is recommended to pay special attention to the following:
- It is recommended that patients undergoing MRI scans are thoroughly monitored for perceived temperature and/or pain sensations.
 - Patients with impaired thermoregulation or temperature sensation should be excluded from MRI scanning procedures.
 - In general, it is recommended that a low field strength MRI system is used when conductive implants are present. The specific absorption rate (SAR) used should be reduced as much as possible.
 - The use of the ventilation system can further help to reduce the temperature rise in the body.

The implants are made of materials that are not magnetic. For this reason, the systems have not been tested for heating, migration, or imaging artifacts in the magnetic resonance environment. The safety of these systems in the magnetic resonance environment is unknown.

Note that MRI techniques in the form of dedicated spin-echo sequences such as VAT, SEMAC (SEMAC-VAT - 2D), MSVAT-SPACE (3D), and MAVRIC should be used to optimize image quality and minimize artifacts

Single-Use Product

The single-use implants are supplied sterile and must not be reprocessed.

The medical devices are gamma sterilized and comply with EN ISO111347 (parts 1 and 2) and EN ISO11607 (parts 1 and 2) standards.

Precautions

- Confirm functionality of instruments
- Select the implants according to the weight, the patient's activity level and the bone fracture to be restored.
- It is recommended to monitor the patients to prevent infections and pain.
- Note that the forces to be transmitted by the implants are kept low by selecting suitable biomechanics. Otherwise, overloading may cause the screws to tear out/break.
- Excessive loading of the implants by the patient's body weight and activity level should be avoided.
- Extreme deformation of the implants must be avoided. Careful bending of the plates, performed with appropriate care, will not damage the implants.
- Repeated deformation must be avoided at all costs to prevent material fatigue and/or fracture of the implant.
- Reuse of implants is absolutely forbidden (single use).
- We strongly advise informing the patient about the advantages and disadvantages of implants.

Not following these precautions can have serious consequences for the healing process.

Possible negative side effects

Risks, side effects, and adverse events may occur with all major surgical procedures

These are:

- Problems resulting from anesthesia and patient positioning (e.g. nausea, vomiting, dental injuries, neurological impairments)

- Thrombosis
- Embolism
- Infection
- Nerve and/or tooth root damage
- Injuries to blood vessels
- Bleeding
- Damage to soft tissues incl. swelling
- Scar formation
- Functional impairment of the musculoskeletal system
- Pain
- Discomfort or abnormal sensation due to the presence of the device
- Allergy
- Side effects associated with the implant (poor or delayed fusion of the bone)

Patient Information

The user should record and keep all information provided to the patient. It should be checked before use whether the patient tolerates the material to be implanted.

Handling

The CRANIOBRIDGE™ implant consists of plates and / or screws.

The implants described in these instructions for use may only be used (implanted) by surgeons with the appropriate experience.

The implants may only be used if the blister packaging is undamaged, otherwise sterility cannot be guaranteed.

After removing the implant carrier from the blister packaging, a visual inspection must be performed before implantation to ensure that all implants are in flawless condition.

In case of non-observance of these instructions or in case of demonstrable violation of the instructions adeor medical AG accepts no responsibility for claims.

Removal of the implants

Normally the implants are removed after successful fracture healing. Dispose the removed implants according to local and national regulations.

Storage

Store in a dry place, protected from dust, without external force, without major temperature fluctuations and away from any aggressive media.

Disposal

Unused implants must be properly disposed in accordance with local and national regulations. Under no circumstances may these be reprocessed and implanted / used at another time.

Shelf life

The shelf life of the product is indicated by the expiry date on the packaging.

Symbols

	Order No.
	Lot number
	Manufacturer
	Date of manufacture
	Sterilized by radiation
	Do not reuse / Single use only

 Use by – symbol is accompanied by a date

 Use only with undamaged packaging

 Follow Instructions for Use

 CE-marking including Notified body code

 Medical Device

Warranty

The products of adeor medical AG are developed and manufactured with greatest care. Tests and inspections by highly qualified personnel guarantee safe and perfect functioning. Please note that warranty claims are only valid if all instructions in the enclosed instructions for use are followed.

adeor medical AG is not liable for damage caused by improper handling or repair by third parties not authorized by adeor.

The use of the instruments is the responsibility of the user. Any liability on the part of adeor medical AG is excluded for any resulting damage.

Warranty claims must be submitted to adeor or the service partner, enclosing the proof of purchase.

The provision of warranty service does not extend the warranty period or any other guarantee period.

adeor medical AG
Martinshof 5
83626 Valley
Deutschland

Tel. +49 8024 47717 0
Fax +49 8024 47717 99

eMail: office@adeor.com
http://www.adeor.com

Version 2022-07-DE_EN

adeor medical AG