IMPLANT PASSPORT

IMPLANTATAUSWEIS

PATIENT DATA | PATIENTENDATEN

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**This person received a medical implant.
This may interfere with metal detection devices.**

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MRI SAFETY INFORMATION

adeor medical AG implants are made of titanium and titanium alloys. This type of non-ferromagnetic material poses no risk to the patient during magnetic resonance imaging (MRI) examinations.

Various studies and test reports have shown that these materials can be safely used in an MRI environment according to ISO 5832-2 and ISO 5832-3 and can be considered „MRI Conditional“ according to ASTM F2503. MRI Conditional is defined as an item that has been shown to pose no known hazards in a specific MRI environment with specified conditions of use (ASTM F2503, 3.1.9). Precautions: The above statement is based on non-clinical testing. The actual patient temperature rise depends on a variety of factors beyond SAR and timing of radiofrequency application.

Therefore, it is recommended to pay special attention to the following: It is recommended that patients undergoing MRI scans are thoroughly monitored for perceived temperature and/or pain sensations. Patients with impaired thermoregulation or temperature sensation should be excluded from MRI scanning procedures. In general, it is recommended that a low field strength MRI system is used when conductive implants are present. The specific absorption rate (SAR) used should be reduced as much as possible. The use of the ventilation system can further help to reduce the temperature rise in the body. The implants are made of materials that are not magnetic. For this reason, the systems have not been tested for heating, migration, or imaging artifacts in the magnetic resonance environment. The safety of these systems in the magnetic resonance environment is unknown. Note that MRI techniques in the form of dedicated spin-echo sequences such as VAT, SEMAC (SEMAC-VAT - 2D), MSVAT-SPACE (3D), and MAVRIC should be used to optimize image quality and minimize artifacts.

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