
Bipolar Forceps Reusable

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EN Instructions for use

US Instructions for use

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
INSTRUCTIONS FOR USE BIPOLAR FORCEPS

English

REF

Product group:
HF Forceps

adeor medical AG
Martinshof 5
83626 Valley
Germany

Attention:

The instruments are supplied non-sterile and must be cleaned thoroughly and, if required, disinfected and sterilized before the first and every subsequent use.

Warning:

Legal regulations stipulate that the instruments may only be acquired by or on behalf of a physician.

Warning:

Please read the information in this leaflet carefully. Improper use and care or use for unintended purposes can lead to premature wear or risks for patients and users.

Intended use:

Bipolar forceps made by adeor medical AG are intended for use in all areas of open surgery. They are designed to grasp, manipulate and coagulate tissue. They are connected through a suitable bipolar cable with the bipolar output of a HF generator and may be used only with bipolar coagulation current. adeor bipolar forceps must be operated with the following parameters: Frequency range between 300 kHz and 1,000 kHz; maximum generator operating voltage 600 V_p.

The surgeon must determine the appropriate type of treatment for each individual case in cooperation with the internist and the anaesthetist. Typical HF surgical interventions include

- Dermatology: hair and wart removal
- ENT: polyp removal, bipolar hemostasis
- Gynaecology: electro-conization of the portio, sterilization by bipolar coagulation of the fallopian tube
- Urology: transurethral resection of the prostate, bladder tumour or bladder wall papillomas; cutting and coagulation in transabdominal surgery
- Laparoscopy and cystoscopy: Removal of polyps
- General surgery: primary sections, separating tissue, plastic surgery, hemostasis, coagulation
- Neurosurgery: Bipolar coagulations

Contraindications:

The bipolar forceps should not be used if the physician determines that the risk to the patient outweighs the benefits.

Bipolar forceps made by adeor medical AG should not be used for tubal sterilization or tube coagulation for sterilization. Do not use on patients with cardiac pacemakers or other active implants.

No or limited HF surgical interventions should take place:

- for patients with implantable cardiac pacemakers or pacemaker electrodes
- in areas where flammable or explosive materials are present, for example in the gastro-intestinal tract (fire and explosion hazard)
- in acute emergency situations
- in the case of severe coagulation disorders
- in the case of severely impaired pulmonary or cardiovascular functions

Attention:

Instruments for electrosurgery may only be used by persons who are specially trained or instructed.

Safety instructions:

- Checking, using and handling electrosurgical instruments is the responsibility of the user.
- Before the first and every other application, each instrument must be cleaned thoroughly, disinfected and sterilized.
- Before each use, the instrument must be inspected for wear and visible damage, such as cracks, fractures or damage to the insulation, in particular areas such as the blade, tip, groove, connector, lock and moving parts. Check the insulation thoroughly.
- Do not use damaged instruments.
- Do not use the instruments in the presence of flammable or explosive substances.
- Do not place the instrument on top of the patient.
- Perform coagulation only if the contact surfaces are in the visible area. Do not touch any other metal instruments during the procedure.

Reprocessing / Number of processing cycles / Service life:

The product life span and number of treatment cycles depend on how carefully the product is handled and whether the product processing instructions are followed. The service life may be reduced by a non-specifiable timeframe depending on which type of applications the product is used for or due to repeated processing. After a maximum of 100 applications, the forceps need to be checked for functionality as well as for the correct state of insulation and connector. If the forceps are damaged, especially at the insulation or the connector, you must dispose of the product or send it to the manufacturer for repair.

Preparation and transport:

The instruments must be stored and transported to the processing location in a closed container to avoid damaging the instruments or contaminating the environment.

Processing**Warning:**

A manual processing procedure was not validated as the product is not suitable on the basis of its design (grooves for grasping and difficult-to-access places at the base).

Automatic cleaning:**Note**

The reprocessing and disinfection procedures are validated according to current valid standards.

Mechanical reprocessing and disinfection

Cleaning and disinfection procedure were validated on a WD BHT INNOVA® M3 device under the following program:

1	Cleaning phase	1 min 30°C / 86°F
2	Cleaning phase	6 min 55°C / 131 °F
3	Rinse	1 min demin. water
4	Thermic disinfection	5 min 90°C / 194°F demin. water

Cleaner:

0.5% neodisher Mediclean Forte (Dr. Weigert).

Check visible surfaces for residues after mechanical cleaning/disinfecting. Repeat the cleaning procedure if necessary.

Drying

In accordance with the program for the washer and disinfectant.

Sterilization / Autoclaving**Packaging**

The product should be single wrapped in sterile pouches (e.g., paper-laminate ASSURE Plus sterilization pouches – foil accord. ISO 11607)

Sterilization process

Steam sterilization using fractionated pre-vacuum procedure. The maximum load for the steriliser may not be exceeded.

Pre-vacuum method: 4 min, 132 °C / 272 °F
 Drying time: the autoclave calculates the time of the drying process

adeor recommends a drying time of 20 min.

The pouch/forceps have to be visibly dry at the end of the sterilization process.

Note:
 The sterilization procedure was validated according to current standards.

Attention:

We do not take any responsibility for reusing tools if the instruments are used for patients with Creutzfeldt-Jakob-, Scrapie- or BSE-infection.

Non-stick Bipolar forceps:

The polished tips are made of silver (with gold plating) and may tarnish which does not affect the function of the instruments. You can restore the original shine by rubbing the non-stick bipolar forceps with a silver polishing cloth.

Functional testing:

The instruments must be checked for functionality before each use. Visually inspect the instruments for cleanliness and carry out a function test. If necessary, repeat the reprocessing process until the instrument looks clean. Surface damage, such as scratches, cracks, nicks, notches etc., or bent parts mean that the instruments need to be repaired and may not be used. Do not use damaged instruments.

Warning:

Malfunction, instrument failure, damaged insulation

Action required by user:

Perform a visual inspection/functional testing. It is the responsibility of the user to ensure

that the instruments are used in line with the intended purpose.

Warning:

Non-sterile residue on the forceps

Action required by user:

Clean before first use and sterilize according to the processing instructions.

Warning:

Electrical loading capacity of the electrodes exceeded

Action required by user:

Follow manufacturer's instructions for electrical loading; see "Intended use".

Warning:

Disinfectant has not evaporated completely

Action required by user:

Follow the instructions for thermal disinfection.

Warning:

Instruments were not cleaned properly

Action required by user:

Clean the instruments according to the processing instructions.

Warning:

Moisture penetration, corrosion

Action required by user:

Follow storage and packaging instructions.

Repairs/service delivery for repairs at adeor medical AG:

Do not perform repairs. Service and repair work must be carried out by trained and qualified personnel only. Please contact the manufacturer for service and repair queries.

Attention: Before you send defective products back for repair, these products must have undergone the entire reprocessing process as described in the processing instructions. Please enclose a statement or

evidence of cleaning, disinfection and sterilization with your return.

Handling:

Surgical instruments should always be treated with the utmost care during transport, cleaning, care, sterilization and storage, in particular pointed tips and other sensitive areas.

Generator connections:

Forceps made by adeor medical AG are approved for use with the following generators:

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Type	Item no.	max. bipolar coagulation power
VIO 300 D	10140-100	≤ 200 Watt
VIO 200 D	10140-200	≤ 120 Watt
VIO 300 S	10140-300	≤ 200 Watt
VIO 200 S	10140-400	≤ 120 Watt
VIO 100 C	10140-500	≤ 80 Watt
VIO 50 C	10140-550	≤ 50 Watt
ICC 350	ICC 350	≤ 150 Watt
ICC 300	ICC 300	≤ 150 Watt
ICC 200	ICC 200	≤ 150 Watt
ICC 80	ICC 80	≤ 50 Watt
ICC 50	ICC 50	≤ 50 Watt
ACC 451	ICC 451	≤ 120 Watt
ACC 450	ICC 450	≤ 120 Watt
ACC 430	ICC 430	≤ 120 Watt

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Type	Item no.	max. bipolar coagulation power
BM-780 II	360080-01	≤ 70 Watt

KARL STORZ GmbH & Co. KG

Type	Item no.	max. bipolar coagulation power
AUTOCON® II 200	205322 20	≤ 120 Watt
AUTOCON® II 400	205352 20	≤ 120 Watt

COVIDIEN

Type	Item no.	max. bipolar coagulation power
Force FX™	Force FX™	≤ 70 Watt
ForceTriad™	ForceTriad™	≤ 95 Watt
Force EZ™	Force EZ™	≤ 70 Watt
SurgiStat™	SurgiStat™	≤ 30 Watt

BOWA-electronic GmbH & Co. KG

Type	Item no.	max. bipolar coagulation power
ARC 400	900-400	≤ 350 Watt
ARC 350	900-351	≤ 350 Watt
ARC 303	900-303	≤ 120 Watt
ARC 250	900-250	≤ 120 Watt
ARC 100	900-100	≤ 100 Watt

Olympus Surgical Technologies Europe

Type	Item no.	max. bipolar coagulation power
ESG-100, 220-240 V	WB991036	≤ 120 Watt
ESG-100, 100-120 V	WB991046	≤ 120 Watt
ESG-400	WB91051W	≤ 320 Watt

Söring GmbH

Type	Item no.	max. bipolar coagulation power
MBC 600	MBC 600	≤ 100 Watt
MBC 601	MBC 601	≤ 100 Watt
MBC 601 UAM	MBC 601 UAM	≤ 100 Watt
MBC 200	MBC 200	≤ 140 Watt
BCC 140	BCC 140	≤ 140 Watt

LAMIDEY NOURY MEDICAL

Type	Item no.	max. bipolar coagulation power
Optima 4	Optima 4	≤ 370 Watt
Optima 3	Optima 3	≤ 370 Watt
Optima 2	Optima 2	≤ 200 Watt
MC 2	MC 2	≤ 70 Watt
MC 3	MC 3	≤ 100 Watt
MC 4	MC 4	≤ 100 Watt

Integra LifeSciences Corporation

Type	Item no.	max. bipolar coagulation power
Elektrotom® 621	Elektrotom® 621	≤ 95 Watt
Elektrotom® 630	Elektrotom® 630	≤ 95 Watt

Gebrüder Martin GmbH & Co. KG

Type	Item no.	max. bipolar coagulation power
Minicutter	80-008-03-04	≤ 70 Watt
ME 102	80-010-02-04	≤ 80 Watt
maxium®, m-Version	80-042-00-04	≤ 320 Watt
maxium®, i- Version	80-042-02-04	≤ 320 Watt
maxium®, e- Version	80-042-04-04	≤ 320 Watt
ME MB 3, m-Version, 220-240 V	80-040-11-04	≤ 100 Watt
ME MB 3, m-Version, 100-127 V	80-040-11-10	≤ 100 Watt
ME MB 3, i- Version, 100-127 V	80-040-12-04	≤ 100 Watt
ME MB 3, i- Version, 100-127 V	80-040-12-10	≤ 100 Watt
ME 411	80-041-01	≤ 80 Watt

Aesculap AG / B. Braun Melsungen AG

Type	Item no.	max. bipolar coagulation power
GN300	GN 300	≤ 80 Watt
GN640	GN640	≤ 100 Watt
GN060	GN060	≤ 50 Watt

Suitable connection cables:

adeor bipolar forceps can be used with the following bipolar connection cables for the above generators:

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Cable connection, forceps
European flat plug
2-pin plug (US male connector)

These cables may be connected only to the bipolar output of electrosurgical devices. Please follow the generator instructions on connecting these cables before using the bipolar forceps. HF generators pose a number of risks, such as incorrect operation, unintended high-frequency burns, ignition of flammable liquids and gases (risk of explosion).

Packing/Storage:

Some parts of these instruments can be easily damaged. They should therefore be stored individually in their packaging or in a protective container with individual compartments.

Take special care that there are no chemicals nearby to avoid corrosion. The instruments must be packaged in accordance with ISO 11607 for sterilized devices. The sterilized instruments should be stored in a dry, clean and dust-free environment at moderate temperatures of 5°C to 40°C.

Warranty:

adeor medical AG provides customers with tested, fault-free products. All our products are designed and manufactured to meet highest quality standards. However, should problems occur, please contact our customer service team. Liability is excluded for products that were modified from the original product, used for other than the intended purposes, handled or used improperly. adeor medical AG does not accept liability for incidental or consequential damages. This warranty is null and void if repairs are carried out by companies that are not authorized for repair by adeor medical AG.

Proper disposal of bipolar forceps:

If the bipolar forceps can no longer be used due to wear or damage, they must be disposed of properly. This means that the forceps need to be disassembled, impurities removed and the instruments must be sterilized again before disposal (see processing instructions).

Symbols and their meaning:

Follow operating instructions

Protect the product against moisture and condensation

Protect the product against direct sunlight

The product is delivered non-sterile

Caution: Please follow instructions

Manufacturer / Date of manufacture

Order number

Lot number

Prescription device only (USA)

CE including Notified Body code


Medical Device

INSTRUCTIONS FOR USE BIPOLAR FORCEPS

US

REF **Product group:**
HF Forceps

adeor medical AG
Martinshof 5
83626 Valley
Germany

Attention:

The instruments are supplied non-sterile and must be cleaned thoroughly and, if required, disinfected and sterilized before the first and every subsequent use.

Warning:

Legal regulations stipulate that the instruments may only be acquired by or on behalf of a physician.

Warning:

Please read the information in this leaflet carefully. Improper use and care or use for unintended purposes can lead to premature wear or risks for patients and users.

Intended use:

The adeor medical Non-stick Bipolar Forceps are intended for use by a physician familiar with electrosurgery in bipolar coagulation for general open surgery where coagulation of soft tissue is needed. adeor bipolar forceps must be operated with the following parameters: Frequency range between 300 kHz and 1,000 kHz; maximum generator operating voltage 600 V_p.

The adeor medical Non-stick Bipolar Forceps have not been shown to be effective for tubal sterilization or tubal coagulation for

sterilization procedures and should not be used for these procedures.

The types of surgery intended include:

- ENT
- Gynecology (except for use in female sterilization)
- Urology
- General surgery
- Neurosurgery
- Laryngeal Surgery
- Orthopedic Surgery
- Thoracic Surgery

Contraindications:

The bipolar forceps should not be used if the physician determines that the risk to the patient outweighs the benefits.

Bipolar forceps made by adeor medical AG should not be used for tubal sterilization or tube coagulation for sterilization. Do not use on patients with cardiac pacemakers or other active implants.

No or limited HF surgical interventions should take place:

- for patients with implantable cardiac pacemakers or pacemaker electrodes
- in areas where flammable or explosive materials are present, for example in the gastro-intestinal tract (fire and explosion hazard)
- in acute emergency situations
- in the case of severe coagulation disorders
- in the case of severely impaired pulmonary or cardiovascular functions

Attention:

Instruments for electrosurgery may only be used by persons who are specially trained or instructed.

Safety instructions:

- Checking, using and handling electrosurgical instruments is the responsibility of the user.
- Before the first and every other application, each instrument must be cleaned thoroughly, disinfected and sterilized.

- Before each use, the instrument must be inspected for wear and visible damage, such as cracks, fractures or damage to the insulation, in particular areas such as the blade, tip, groove, connector, lock and moving parts. Check the insulation thoroughly.
- Do not use damaged instruments.
- Do not use the instruments in the presence of flammable or explosive substances.
- Do not place the instrument on top of the patient.
- Perform coagulation only if the contact surfaces are in the visible area. Do not touch any other metal instruments during the procedure.

Reprocessing / Number of processing cycles / Service life:

The product life span and number of treatment cycles depend on how carefully the product is handled and whether the product processing instructions are followed. The service life may be reduced by a non-specified timeframe depending on which type of applications the product is used for or due to repeated processing. After a maximum of 100 applications, the forceps need to be checked for functionality as well as for the correct state of insulation and connector. If the forceps are damaged, especially at the insulation or the connector, you must dispose of the product or send it to the manufacturer for repair.

Preparation and transport:

The instruments must be stored and transported to the processing location in a closed container to avoid damaging the instruments or contaminating the environment.

Processing

Warning:

A manual processing procedure was not validated as the product is not suitable on the basis of its design (grooves for grasping and difficult-to-access places at the base).

Automatic cleaning:

Note

The reprocessing and disinfection procedures are validated according to current valid standards.

Mechanical reprocessing and disinfection

Cleaning and disinfection procedure were validated on a WD BHT INNOVA® M3 device under the following program:

1	Cleaning phase	1 min 30°C / 86°F
2	Cleaning phase	6 min 55°C / 131 °F
3	Rinse	1 min demin. water
4	Thermic disinfection	5 min 90°C / 194°F demin. water

Cleaner:

0.5% of a cleaning detergent containing less than 5% nonionic and anionic surfactants, pH value 10.4 – 10.8 (e.g., Neodisher Mediclean Forte Dr. Weigert).

Check visible surfaces for residues after mechanical cleaning/disinfecting. Repeat the cleaning procedure when necessary.

Drying

In accordance with the program for the washer and disinfectant.

Sterilization / Autoclaving

Packaging

The product should be single wrapped in sterile pouches. (e.g. paper-laminate ASSURE Plus sterilization pouches – foil accord. ISO 11607)

Sterilization process

Steam sterilization using fractionated pre-vacuum procedure (sterilization device compliant to EN 285 in its current revision and validated according to EN ISO 17665-1 in its current revision).

Gravity method: 15 min, 132 °C / 272 °F

Drying time: the autoclave calculates the time of the drying process

Pre-vacuum method: 4 min, 132 °C / 272 °F
Drying time: the autoclave calculates the time of the drying process

adeor recommends a drying time of 20 min.

The pouch/forceps have to be visibly dry at the end of the sterilization process.

Note:
The sterilization procedure was validated according to current standards.

Additional Information

The maximum load for the steriliser may not be exceeded.

Attention:

We do not take any responsibility for reusing tools if the instruments are used for patients with Creutzfeldt-Jakob-, Scrapie- or BSE-infection.

Non-stick Bipolar forceps:

The polished tips are made of silver (with gold plating) and may tarnish which does not affect the function of the instruments. You can restore the original shine by rubbing the non-stick bipolar forceps with a silver polishing cloth.

Functional testing:

The instruments must be checked for functionality before each use. Visually inspect the instruments for cleanliness and carry out a function test. If necessary, repeat the reprocessing process until the instrument looks clean. Surface damage, such as scratches, cracks, nicks, notches etc., or bent parts mean that the instruments need to be repaired and may not be used. Do not use damaged instruments.

Warning:

Malfunction, instrument failure, damaged insulation

Action required by user:

Perform a visual inspection/functional testing. It is the responsibility of the user to ensure that the instruments are used in line with the intended purpose.

Warning:

Non-sterile residue on the forceps

Action required by user:

Clean before first use and sterilize according to the processing instructions.

Warning:

Electrical loading capacity of the electrodes exceeded

Action required by user:

Follow manufacturer's instructions for electrical loading; see "Intended use".

Warning:

Disinfectant has not evaporated completely

Action required by user:

Follow the instructions for thermal disinfection.

Warning:

Instruments were not cleaned properly

Action required by user:

Clean the instruments according to the processing instructions.

Warning:

Moisture penetration, corrosion

Action required by user:

Follow storage and packaging instructions.

Repairs/service delivery for repairs at adeor medical AG:

Do not perform repairs. Service and repair work must be carried out by trained and qualified personnel only. Please contact the manufacturer for service and repair queries.

Attention: Before you send defective products back for repair, these products must have undergone the entire reprocessing process as described in the processing

instructions. Please enclose a statement or evidence of cleaning, disinfection and sterilization with your return.

Handling:

Surgical instruments should always be treated with the utmost care during transport, cleaning, care, sterilization and storage, in particular pointed tips and other sensitive areas.

Suitable connection cables:

adeor bipolar forceps can be used with bipolar connection cables that are suitable for use with generators that meet the required operating parameters.

Cable connection, forceps
European flat plug
2-pin plug (US male connector)

These cables may be connected only to the bipolar output of electrosurgical devices. Please follow the generator instructions on connecting these cables before using the bipolar forceps. HF generators pose a number of risks, such as incorrect operation, unintended high-frequency burns, ignition of flammable liquids and gases (risk of explosion).

Packing/Storage:

Some parts of these instruments can be easily damaged. They should therefore be stored individually in their packaging or in a protective container with individual compartments.

Take special care that there are no chemicals nearby to avoid corrosion. The instruments must be packaged in accordance with ISO 11607 for sterilized devices. The sterilized instruments should be stored in a dry, clean and dust-free environment at moderate temperatures of 5°C to 40°C.

Warranty:

adeor medical AG provides customers with tested, fault-free products. All our products are designed and manufactured to meet

highest quality standards. However, should problems occur, please contact our customer service team. Liability is excluded for products that were modified from the original product, used for other than the intended purposes, handled or used improperly, adeor medical AG does not accept liability for incidental or consequential damages. This warranty is null and void if repairs are carried out by companies that are not authorized for repair by adeor medical AG.

Proper disposal of bipolar forceps:

If the bipolar forceps can no longer be used due to wear or damage, they must be disposed of properly. This means that the forceps need to be disassembled, impurities removed and the instruments must be sterilized again before disposal (see processing instructions).

Symbols and their meaning:

Follow operating instructions

Protect the product against moisture and condensation

Protect the product against direct sunlight

The product is delivered non-sterile

Caution: Please follow instructions

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